PDS No. 227250	PRODUCT DATA SHEET	Page 1 of 1
Davisian 00	OncoQuick [®]	6
Revision 09	Greiner Item-No. 227250	greiner
Valid for Item-No.:	227250 (aseptic)	

1.	Description / Specification			
1.1	Description	Consists of a tube with conical bottom and graduation, screw cap, porous barrier 227250: pre-filled with separation medium, aseptic		
1.2	Dimensions	See Customer Drawing		
1.3	Volume	Total volume: 50 ml Sample volume: 15 – 30 ml		
1.4	Material / Resin	Tube: PP (Polypropylene), free of heavy metal Porous barrier: HDPE (High Density Polyethylene), free of heavy metal Screw cap: HDPE (High Density Polyethylene), free of heavy metal Separation Medium: - contains trypan blue - not free of endotoxin		
1.5	Colour	Tube: translucent; print of graduation: blue Porous barrier: white Screw cap: blue Separation Medium: blue		
1.6	Sterilization	No (aseptic)		
1.7	Quality Control	Raw Material-Control: physical testing Product-Control: testing of attributive and variable characteristics in accordance with the valid specification		
1.8	Intended Use	General laboratory product for the enrichment of disseminated circulating tumor cells from whole blood to be used by qualified personnel in a laboratory environment.		
1.9	Other Information	- for single use only - protect from direct light - Instruction manual: see package or homepage: https://www.gbo.com/fileadmin/imported from old/Downloads/Product Flyers/Product Flyers BioScience/999999 UserGuide OncoQuick E.pdf		

2.	Features	
2.1	Basic features	-
2.2	Temperature range	+4°C to +25° (storage temperature)
		+18°C to +25° (working temperature)
2.3	Autoclavability	No
2.4	Centrifugation, max. RCF	1.600 x g: swinging-bucket rotor
2.5	Chemical Resistance	See homepage:
		https://www.gbo.com/en_INT/know-how-services/download-center.html
2.6	Shelf life	Depending on the indication of the OncoQuick® separation medium –
		manufacturer, maximum 2 years
2.7	Other Information	-

3.	Packaging	227250
3.1	Pieces / Rack / Bag	10
3.2	Pieces / Box	10
3.3	Lot-No.	E YY MM XXX (manufacturing facility, year, month, consecutive SAP-No.)
3.4	Other Information	Instruction manual / box

4.	Other Information	
	Research use only. Not for diagnostics.	

Data Sheet subject to change without notice!

Prior Issue	Drawn	Approved	Released	CONFIDENTIAL: Information contained in this
Revision	Date	Date	Date	document or drawing is confidential and proprietary to Greiner Bio-One GmbH. This
08	23 January 2025	24 January 2025	24 January 2025	document may not be reproduced for any
Date	Name	Name	Name	reason without written permission from Greiner Bio-One GmbH. All rights of design, invention,
24.06.2024	S. Kaelberer	Dr. T. Schreiber	Dr. CK. Chai	and copyright are reserved.